



Easy
Grip



Easy
Application



Leak
Proof

An easy and convenient way to irrigate patient wound sites. The device consists of a hard plastic syringe with a soft rubber bulb that is squeezed to commence irrigation. Available in sterile and non-sterile versions, suitable for a wide variety of clinical settings.

EASY TO READ MEASUREMENT

60cc barrel with 5cc graduations and half ounce graduation from ½ to 2oz

EASY GRIP BULB

Controls drawing or expelling with use of only one hand

LEAK PROOF

For maximum control in your clinical setting



Suitable for:

Theatre

First aid

Reorder	Size	Sterility	Pcs/Unit	Units/Ctn
BSY001	60ml	Non-Sterile	50	1
BSY100	60ml	Sterile	50	1