



Waterproof  
Backing



Moisture Vapour  
Permeable



Contours  
to Body

Osmocel Hydroporous Foam is an ideal dressing to manage wounds with different needs and moisture levels.

The fine pores retain moisture at the wound surface, while excess exudate travels into upper pores of the foam to be eliminated as moisture vapour, increasing dressing wear time.

## FINE POROUS INTERFACE

Absorbs excess exudate while retaining moisture at wound surface

## THERMAL INSULATION

Creates the ideal temperature to encourage cell division and migration

## LOW DENSITY

Decreasing the risk of profiling and indentations from the dressing

## EXTENSIVE RANGE

Sacral and vari-site options for accurate fit



## Suitable for:

Primary and/or secondary dressing for the management of wounds

Arterial, diabetic and venous leg ulcers

Pressure injuries, partial thickness burns and surgical incisions

Reorder	Size	Pcs/Unit	Units/Ctn
<b>STANDARD</b>			
<b>AWF050</b>	5cm x 5cm	10	10
<b>AWF051</b>	6cm x 8.5cm	5	10
<b>AWF052</b>	10cm x 10cm	10	10
<b>AWF053</b>	10cm x 20cm	5	10
<b>AWF054</b>	15cm x 15cm	5	10
<b>AWF055</b>	20cm x 20cm	10	10

Reorder	Size	Pcs/Unit	Units/Ctn
<b>VARI-SITE</b>			
<b>AWF180</b>	12cm x 19cm	5	10
<b>SACRAL</b>			
<b>AWF150</b>	18cm x 18cm	5	10
<b>AWF151</b>	23cm x 23cm	5	10
<b>SACRAL BUTTERFLY</b>			
<b>AWF160</b>	18cm x 18cm	5	10
<b>AWF161</b>	23cm x 23cm	5	10