

ThreadX® X-Ray Detectable Peanut Swab in Tray



Highly Absorbent



100% Cotton



X-Ray Detectable

ThreadX Peanut Swabs are constructed from 100% cotton and are highly absorbent with a black x-ray detectable thread heat bonded to the gauze. Supplied in a plastic tray of 5 with an adhesive bottom for ease of use.

ADHESIVE TRAY


Prevents tray from slippage or movement

X-RAY DETECTABLE

X-ray detectable thread for patient safety

MOULDED TRAY

Provides quick and visible countability



Suitable for:

Open surgical procedures involving delicate anatomy

Absorption of exudate during neurosurgery

Reorder	Pcs	Sterility	Pcs/Unit	Units/Ctn
XPE001	5pcs	Sterile	50	1

